



Millennia **Command**TM


Simplify workflow and improve performance with the new **Millennia Command**TM slide: the slide that checks all the boxes.

Strong adhesion meets the rigorous demands of harsh IHC processing, while low background staining makes Millennia Command an excellent choice for special stains and H&E; a true all-purpose slide.

Millennia Command

Hybrid slide suitable for IHC, H&E and special staining. Optimized coating for strong adhesion to reduce tissue loss and provide consistency of staining, with minimal background staining during all types of staining. Hydrophilic.

Millennia Command, white frosted, 100/box	MCOMM/W/90
Millennia Command, pink frosted, 100/box	MCOMM/P/90
Millennia Command, green frosted, 100/box	MCOMM/G/90
Millennia Command, blue frosted, 100/box	MCOMM/B/90
Millennia Command, yellow frosted, 100/box	MCOMM/Y/90
Millennia Command, red frosted, 100/box	MCOMM/R/90
Millennia Command, tan frosted, 100/box	MCOMM/T/90
Millennia Command, lilac frosted, 100/box	MCOMM/L/90



Why switch to Millennia Command?

High adhesion produces excellent IHC results

Low background staining on all applications

Validated on Leica and Ventana platforms

Compatible with all thermal transfer printers

- StatLab ESPO
- SlideMate and SlideMateAS
- Primera
- Leica HistoCore PERMA
- TissueTek SmartWrite

Side by side studies reveal equivalency compared to popular brand in:

- ☞ Stain Quality
- ☞ Stain Uniformity
- ☞ Tissue Integrity

General Data
HEALTHCARE