

# FWS™ & TFB™

## Smart Waterbaths



**FWS-S**  
Low-profile deep dish  
waterbath



**FWS-L**  
Waterbath shown with front-mounted  
slide drying rack (SDR-FL)



**TFB-S**  
Low-profile deep dish  
waterbath



**TFB-L**  
Waterbath shown with side-mounted  
slide drying rack (SDR-S)

### ■ Choose the Size that Best Suits Your Needs

The small footprint FWS-S and TFB-S are ideal for crowded workspaces; the deep dish design allows for efficient ribbon floating and pickup. The FWS-L and TFB-L have an 8" x 11" removable glass dish. The FWS models feature built-in heated slide drying areas and a 2" x 2" HISTOorientor for removing wrinkles. All models have heat-resistant plastic housing and a low profile for user comfort.

### ■ Microprocessor-Controlled Temperature Regulation

Temperature is precisely controlled to prevent tissue imperfections such as cracking and layer separation that can result from temperature fluctuation. *All heated areas (dishes, drying areas, and HISTOorientators) can be individually temperature-controlled according to user preference.*

### ■ Optimum Lighting with No Glare

The LED light array combined with a high-contrast background allows greater visibility than conventional waterbaths. The light array is strategically positioned on the side of the glass dish (away from the user) so maximum light is shed on ribbons with the least amount of glare.

### ■ Simple to Clean

Glass dishes are easily removed, and all surfaces are easily accessible and chemical resistant to stand up to frequent cleaning.

**Protect every life story™** in your lab.

# Product Specifications:

8"x 11" Model Dimensions: WxLxH(in)(cm)	(15 x 18.75 x 5.25) (38.1 x 47.6 x 13.3) (in) (cm)
8"x 8" Model Dimensions: WxLxH(in)(cm)	(11.75x18.75x5.25) (29.8x47.6x13.3) (in)(cm)
Weight :lbs(kg)	11 lbs; 6kg
Voltage:	100-120 VAC, 50/60Hz, 5A and 220-240VAC, 50/60Hz, 2.4A
Certification:	TUV/FCC

Catalog #	Description
FWS-L	Waterbath, microprocessor controlled w/11" x 8" removable glass bowl w/built-in LED light array, 2" x 8" heated slide dryer and 2"x 2" HISTO <b>Oriente</b> r™, 100-120VAC, 50/60Hz, 5A.
FWS-L-220	Waterbath, microprocessor controlled w/11" x 8" removable glass bowl w/built-in LED light array, 2" x 8" heated slide dryer and 2"x 2" HISTO <b>Oriente</b> r™, 220-240VAC, 50/60Hz, 2.4A.
FWS-S	Waterbath, microprocessor controlled w/8" x 8" x 2 ¼" deep dish removable glass bowl w/built-in LED light array, 2" x 6" heated slide dryer and 2"x 2" HISTO <b>Oriente</b> r™, 100-120VAC, 50/60Hz, 5A.
FWS-S-220	Waterbath, microprocessor controlled w/8" x 8" x 2 ¼" deep dish removable glass bowl w/built-in LED light array, 2" x 6" heated slide dryer and 2"x 2" HISTO <b>Oriente</b> r™, 220-240VAC, 50/60Hz, 2.4A.
TFB-L	Waterbath, microprocessor controlled w/11" x 8" removable glass bowl, built-in clock & LED light array. 100-120VAC, 50/60Hz, 5A.
TFB-L-220	Waterbath, microprocessor controlled w/11" x 8" removable glass bowl, built-in clock & LED light array. 220-240VAC, 50/60Hz, 2.4A.
TFB-S	Waterbath, microprocessor controlled w/8" x 8" x 2 ¼" deep dish removable glass bowl, built-in clock & LED light array. 100-120VAC, 50/60Hz, 5A.
TFB-S-220	Waterbath, microprocessor controlled w/8" x 8" x 2 ¼" deep dish removable glass bowl, built-in clock & LED light array. 220-240VAC, 50/60Hz, 2.4A.

## More Information

844.643.1129

[www.general-data.com/hc](http://www.general-data.com/hc)

*General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.*

©2018 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.  
4043 McMann Road, Cincinnati, OH 45245  
844.643.1129

Email: [talktous@general-data.com](mailto:talktous@general-data.com)  
Web: [www.general-data.com/healthcare](http://www.general-data.com/healthcare)

Service:  
PH: 800.245.5029  
Email: [customersupport@general-data.com](mailto:customersupport@general-data.com)