

Medical Face Masks

Made in the USA



■ Made in the USA

You can trust these medical-quality face masks will live up to your expectations. Even better; no worries about extended lead times or overseas freight issues, as they are made right here in the United States.

■ High Quality 3-Ply Masks

With a three-layer construction, these quality masks have an ASTM International (formerly American Society for Testing and Materials) rating of ASTM Level 3, a BFE (Bacterial Filtration Efficiency) rating of 99.5%, and a Delta P (Breathability) Rating of 3.5 – 3.7. Filters 99% of viruses, bacteria, dust, chemicals, particles, pollen, and smoke.

■ Comfortable and Lightweight, Forms to Face Contours

Latex-free, lightweight and soft, with foam-based ear straps that are easy on the delicate skin around the ears. Simply slip straps over the ears, lightly press along the nose piece to fit, and adjust the mask to cover the nose to the jaw.

■ Designed for Single Use PPE

Our masks are disposable and are for single use only. One size fits most.


Product Specifications:

Dimensions when flat:	(6.90 x 3.75 in.) [17.5 x 9.5 cm]
Color:	Light blue
Packaging	50 masks/box; 40 boxes/carton; total 2000 masks/carton
Material:	Non-woven from polypropylene
Construction:	3 layer (Spunbond-Meltblown-Spunbond 25 GSM (grams per square meter)

FDA 2020 Registration Number:	3016667418
Carton Dimensions:	16 in. (L) x 16 in. (W) x 20 in. (H)
Carton Weight:	20 lbs. net
Carton Weight:	21 lbs. gross
Pallet w/24 cartons total weight:	532 lbs.
Pallet x24 cartons dimensions:	48 in. x 40 in. x 64 in.

PRODUCT TEST RESULTS			
Test Program	Test Standards	Average	Unit
BFE	ASTM 2101	≥ 99.9	%
PFE	ASTM 2299	≥ 99.5	%
Differential Pressure	EN 14683	4.38	mmH ₂ O/min
Flammability	16 CFR Part 1610	IBE (no time)	S
Synthetic Blood	ASTM F 1862	160	mmHg

Catalog #	Description
M2-2MS-9958	Face Masks, 50/box; 40 boxes/carton (Minimum order 5 boxes)



More Information

844.643.1129 www.general-data.com

©2020 General Data Company. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Company, Cincinnati, Ohio, USA.



Smarter Technology. Exceptional Results.

General Data Company
4354 Ferguson, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com

Web: www.general-data.com

Service:

PH: 800.245.5029

Email: customersupport@general-data.com