

Histology Innovation for a NEW Generation

Software

LabCycle

Specimen Identification Solution



**Proven Technology. Superior Reliability.
Protecting Every Life Story In Your Lab.**

Each tissue specimen your laboratory processes has a life story behind it. Your mission is to ensure all your specimens receive optimal processing in order to deliver the best possible care. With ID/Positive cassette markers, slide labelers, barcode scanners and software, laboratories have everything they need to positively and reliably identify, track and manage specimens throughout their entire laboratory process, start to finish. Laboratories trust ID/Positive to eliminate specimen misidentification, streamline laboratory workflow, increase productivity and throughput, and enhance overall patient safety.

LabCycle - Scalable Specimen Identification, Tracking and Reporting Software

LabCycle is a scalable software solution used to identify specimens throughout the histopathology laboratory's processes. In its simplest form, LabCycle provides labeling for specimen requisitions, containers, cassettes and slides. When utilizing the solution to its full potential, LabCycle tracks data at each level of specimen processing, documenting individual interaction with each barcoded specimen along with time stamps and locality. LabCycle's built-in reporting and dashboards display real-time allocation of laboratory resources, while the simple configuration tools allow laboratories to customize the amount of specimen data to capture.

User-Friendly

- Intuitive user screens make it easier to learn the software minimizing the training period
- Customizable prompts allow specific laboratory terms to be utilized
- Role-based access allows users to only view their assigned modules

Scalable

- Laboratories decide how much of the system to implement
- Supports multiple sites
- Modules can be implemented in phases

Specimen Information

- Real-time access to a specimen's current location
- Ability to view a specimen's history
- Tracking of instruments used to process the specimen

Dashboard and Reports

- Easily view the number of specimens processed
- View key performance indicators (KPIs) such as turnaround time (TAT)
- Locate specimens that have not moved throughout the laboratory within an expected timeframe


Hardware Compatibility

- Option to utilize existing modern computer systems (meeting minimum system requirements)
- Native support for existing ID/Positive barcode scanners, slide labelers and laser cassette markers
- Ability to use compatible barcode scanners, slide labelers and cassette printers



LabCycle - Includes Real-Time Reporting of Every Item's Status & Last Known Location

LabCycle is part of General Data's ID/Positive family of barcode specimen identification, tracking and workflow management solutions for laboratories. With ID/Positive cassette markers, slide labelers, barcode scanners and LabCycle, laboratories will have everything they need to positively and reliably identify, track and manage specimens throughout their entire laboratory process, start to finish. Laboratories trust ID/Positive to eliminate specimen misidentification, streamline laboratory workflow, increase productivity and throughput, and enhance overall patient safety.



Benefits of LabCycle

- + System can be interfaced to the Laboratory Information System (LIS), or data can be entered directly into the system
- + Ability to select the modules that would most benefit the laboratory
- + Identifies and tracks specimens
- + One-to-one verification minimizes human error
- + Enhances patient safety
- + Provides the location and history of each specimen


More Information

PH: 844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2014 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
Headquarters: 4354 Ferguson Drive, Cincinnati, OH 45245
TBS Facility: 3014 Croasdaile Drive, Durham, NC 27705

Sales Inquiries:
PH: 844.643.1129
Email: talktous@general-data.com
Web: www.general-data.com/hc

Service:
PH: 800.245.5029
Email: customersupport@general-data.com