

Protect every life story™ in your lab.

LaserTrack™

AB1 Cassette Printer



- **Reliable, High Resolution Text & Barcodes**

The precision laser used in the AB1 cassette printer creates the crisp edges and high resolution, high-contrast alphanumeric text and 2D lines a barcode requires to produce a 99.5% scan rate throughout the tracking process.

- **Easily Removable/Interchangeable Cassette Magazines**

Easily interchangeable magazines hold up to 70 cassettes each.

- **Smooth Operation**

Cassettes load easily into the magazines, and the optical sensor system and simple marking mechanism greatly minimize jamming.

- **Compact Footprint**

At 10"W x 18"D x 25", this unit fits into almost any size work space.



Product Specifications:

Dimensions: WxDxH	(10 x 18 x 25) (25.4 x 45.72 x 63.5) (in) (cm)	Environmental	Normal laboratory temperatures/humidity
Weight: lbs, kg	39 lbs; 19.5 Kg	Connectivity	USB/Ethernet/Serial (to PC)
Voltage	110/220 VAC 50/60 Hz	Barcode Symbologies	Code 128, Code 39, Datamatrix, PDF417 (truncated)
Resolution	200 lpi	Graphics	General scaled/monochrome bitmap only
Speed	5 seconds cassette to cassette (approx. 12 cassettes per minute)	Media Max Load	1 stack of 70 cassettes each

Capable of interfacing with:

- Cerner
- Epic
- Sunquest
- MEDITECH
- Other LIS/LIMS systems

G1-205-CL01 LaserTrack™ AB1 Laser Cassette Printer

Includes:

- Power Cord
- USB Cable
- User Manual
- 1 magazine
- Operating Software

G1-2EE-9844 LaserTrack™ Collection Tray



Cassettes

White Yellow Pink Aqua Blue Orange Gray Green Lilac Peach Tan

Tissue Cassettes

G1-CTP-100(color) (WE) (YW) (PK) (AQ) (BE) (OE) (GY) (GN) (LC) (PH) (TN)

Biopsy Cassettes

G1-CTP-100(color) (WE) (YW) (PK) (AQ) (BE) (OE) (GY) (GN) (LC) (PH) (TN)

Ordering Example: to order Green Tissue Cassettes, the part number would be G1-CTP-100GN

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com