

 **POLYFin**[™]
Paraffin Wax



■ **Dependable, Consistent Formulation**

This unique formulation has ideal handling characteristics; it is less sticky, less tacky, and less static, all important for better handling. No pre-treating is necessary to soften tissue prior to cutting, and since it has a lower melting point (60° C), it will not overcook or harden tissue. Use POLYFin for proper handling and great tissue support.

■ **High Polymer Content**

More polymer means more cutting rigidity, which is important for producing the most precise cutting at the microtomy station. With labs being asked to do more and more quantitative work, precision and efficiency are more important to workflow than ever before. PolyFin provides superior tissue infiltration and texture, allowing for wrinkle-free sections at 4-5 microns with no compression.

For Both Processing and Embedding

- POLYFin Paraffin Wax can be used in both tissue processing and embedding; offering both user convenience and savings.

Protect every life story[™] in your lab.

Product Specifications:

Catalog #	Description
H-PF	Size: 8x1kg (1kg bag, 8 bags/case (17.6 lbs. total))

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com