

Personal ID® AC Color Identification Wristbands



■ Quick & Easy to Use

Error-free patient identification is one of the most important steps in ensuring patient safety. Just type the information, insert the photo, and print. Personal ID barcoded patient wristbands help prevent dangerous medical errors by making patient information immediately available for medication administration, specimen collection, surgical procedures, and any other patient care activity.

■ On-Demand Color Printing

Personal ID AC Color wristbands enable the user to print any color directly on the wristband as the wristband is printed. Only one wristband is needed - a patient is no longer required to wear multiple wristbands for different color-code alerts.

■ Reliable & Durable

Personal ID AC Color wristbands are extremely durable and will stand up to hand sanitizers, soapy & foam washes, exposure and wear of a typical hospital stay.

■ Scans First Time, Every Time

When printed with Epson's TM C3500 printer, they produce high resolution linear and two-dimensional barcodes that are high-contrast and deliver immediate first-time scans.

■ Ultimate in Patient Comfort

Soft, lightweight, and smooth to the patient's skin, with a unique bracelet design in sizes for patients from infant/pediatric to standard adult sizing. They are latex-free and hypo-allergenic, with no skin-irritating agents or chemical additives.



Positive Patient Identification starts here.

Product Specifications:

Adult AC Color		Pediatric AC Color	
Part Number:	L-WB-AIJ29330	Part Number:	L-WB-AIJ29227
Size:	1.5" W x 13.25" L	Size:	1.2" W x 9.12" L
Quantity per Roll:	50 (2" core)	Quantity per Roll:	50 (2" core)
Rolls per Carton:	5	Rolls per Carton:	5
Print Method:	Inkjet	Print Method:	Inkjet
Thickness:	10 Mils	Thickness:	10 Mils

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps facilities enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com