

## Reagents for Processing & Staining



Top quality reagents for reproducible results; providing assurance that you can process and stain tissues of all types. All alcohols are ACS grade; with color-coded labels on each product for error reduction.

- **Formalin, 10% Neutral Buffered Formalin**  
Ready to use phosphate-buffered formalin.
- **Reagent Alcohol 100% (200 proof)**  
Consists of 95 parts Formula 3A specially denatured ACS grade alcohol and 5 parts ACS grade isopropanol. Ideal for virtually all processing and staining procedures in the lab. Ultra-high filtration and quality control ensures consistently high performance.
- **Reagent Alcohol 95%**  
Reagent Alcohol 100% with 5 parts deionized water.
- **Reagent Alcohol 80%**  
Reagent Alcohol 100% with 20 parts deionized water.
- **Reagent Alcohol 70%**  
70% Alcohol blend using an ultrapure type 1 water. Specifically formulated for use with cytology imaging systems, this alcohol blend also performs well in all standard laboratory applications.
- **Xylene**  
Purified and filtered for a consistent, quality product.
- **XydeClear Xylene Substitute**  
Higher permissible exposure limit than xylene. 100% aliphatic hydrocarbon for processing and staining. Virtually no odor.

**Protect every life story™** in your lab.

# Product Specifications:

Catalog #	Description
F-10NB1G	10% Neutral Buffered Formalin, 4 gal/cs.
F-10NB5G	10% Neutral Buffered Formalin, 5 gallon cube
A-100RA	Reagent Alcohol 100% - 200 proof, 4 gal/cs.
A-95RA	Reagent Alcohol 95% - 4 gal/cs.
A-80RA	Reagent Alcohol 80% - 4 gal/cs.
A-70RA	Reagent Alcohol 70% - 4 gal/cs.
C-X1G	Xylene, 4 gal/cs.
C-X3S1G	XydeClear Xylene Substitute, 4 gal/cs.

## More Information

844.643.1129

[www.general-data.com/hc](http://www.general-data.com/hc)

*General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.*

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.  
4043 McMann Road, Cincinnati, OH 45245  
844.643.1129

Email: [talktous@general-data.com](mailto:talktous@general-data.com)  
Web: [www.general-data.com/healthcare](http://www.general-data.com/healthcare)

Service:  
PH: 800.245.5029  
Email: [customersupport@general-data.com](mailto:customersupport@general-data.com)