

SHURDry™ SD-III

Slide Dryer



- **Operator Defined Programs Control Drying Times**
The control panel utilizes seven keys that function as individual timers for multi-user drying times. Keys can be programmed for different times and a distinct alarm will sound when an item is ready to be removed.
- **See-Through Doors Allow Easy Content Monitoring**
Glass doors and sides provide the histologist the visibility needed for examining contents without releasing heat by opening doors. They additionally allow for convenient cleaning.
- **Small Countertop Footprint; Large Internal Capacity**
An adjustable shelf allows for an assortment of sizes of bowls or containers.
- **Keeps Paraffin Liquefied**
The SD-III design makes storage of melted paraffin highly efficient and eliminates the need to re-melt the paraffin.

Protect every life story™ in your lab.

Product Specifications:

External Dimensions: WxDxH (in) [cm]	(22 x 18 x 24.5) [56 x 46 x 62]
Internal Dimensions: WxDxH (in) [cm]	(21x11x14)[52x27x37]
Weight: (lbs) [kg]	(70) [32]
Voltage:	100-120 VAC, 50/60Hz, 6.2A and 220-240VAC, 50/60Hz, 3A

Catalog #	Description
SD-III-120	High-Capacity Slide Dryer/Oven/Incubator; forced air; w/ two shelves (one adjustable); six operator-definable programs for controlling drying times; microprocessor controlled; ambient to 75°C; 9W recessed fluorescent light; four 3,200 rpm fans, 2 speed fan selectability; 700w, 18" heating strip. Specifications: 100-120V, 60Hz, 6.2A.
SD-III-220	High-Capacity Slide Dryer/Oven/Incubator; forced air; w/ two shelves (one adjustable); six operator-definable programs for controlling drying times; microprocessor controlled; ambient to 75°C; 9W recessed fluorescent light; four 3,200 rpm fans, 2 speed fan selectability; 700w, 18" heating strip. Specifications: 220-240V, 50/60Hz, 3A.

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com