

# SHURClean™

Fume Hoods



**SC-36S-X or -F**  
Sloped Front | Width: 36"  
Filters: Xylene or Formalin



**SC-36S-WOME-TV or -RV**  
**SC-48S-WOME-TV or -RV**  
6" Ports: Top or Rear Mounted | Widths: 36" or 48"

## ■ Choose the Style that Best Suits Your Needs

All units are sloped-front fume hoods; feature vibrant LED lighting inside the work area, and utilize hinged sashes for manual grossing and/or staining applications. All units meet operator face velocity requirements of published industry standards with sashes closed.

## ■ Built-In Filtration for Easy Placement

The SC-36S-X/F and SC-DD models have a built-in updraft filter for areas where external HVAC ventilation is not available or possible. This is an ideal solution for adding fume hoods to existing labs without the need for extensive equipment reconfiguration. Filters are available for both xylene and formalin fumes.

## ■ User-Friendly Control Panel

The SC-36S-X/F models feature an LED control panel that enables the user to easily program and control ON/OFF times and fan speeds for optimal efficiencies and fume management.

## ■ Top or Rear Mounted Exhaust Ports

For installations where external ventilation is available, the SC-36S-WOME-TV or -RV, SC-48S-WOME-TV or -RV sloped front models are available with 6" exhaust ports that can be on top or rear mounted.

# Product Specifications:

<b>Active Slope Front Fume Hood</b>	<b>SC-36S-F or -X</b>		
External Dimensions: WxDxH(in)[cm]	(36x31x28) [91x79x71] (in)[cm]		
Internal Dimensions: WxDxH(in)[cm]	(34x19x19) [86x48x48] (in)[cm]		
Front Opening Height: (in)[cm]	<b>Sashes Closed</b> (8) [20]	<b>1 Sash Open</b> (13) [32]	<b>2 Sashes Open</b> (19)[49]
Weight: (lbs)[kg]	(67 lbs) [31 kg]		
Voltage:	100-120 VAC, 50/60Hz, 1.1Amps and 220-240VAC, 50/60Hz, 0.7Amps		
Speed:	Two variable speed 174 cfm blowers	Face velocity: 90 fpm with sash completely down	
Bulb:	22" LED 12VDC, 4 Amps		

<b>Passive Slope Front Fume Hood</b>	<b>SC-36S-WOME-TV or -RV</b>	<b>SC-48S-WOME-TV or -RV</b>
External Dimensions: WxDxH(in)[cm]	(36x22x22) [91x57x55] (in)[cm]	(48x22x22) [122x57x55] (in)[cm]
Internal Dimensions: WxDxH(in)[cm]	(34x20x19) [86x51x49] (in)[cm]	(46x20x19) [117x51x49] (in)[cm]
Weight: (lbs)[kg]	(31 lbs) [14 kg]	

Catalog #	Description
SC-36S-F or -X	Active Sloped-Front Ductless Fume Hood for grossing (-F); Active Sloped-Front Ductless Fume Hood for coverslipping (-X)
SC-36S-WOME-TV or -RV	36" Passive Sloped-Front Ductless Fume Hood for grossing or coverslipping, double-hinged sash
SC-48S-WOME-TV or -RV	48" Passive Sloped-Front Ductless Fume Hood for grossing or coverslipping, double-hinged sash

## More Information

**844.643.1129**  
[www.general-data.com/hc](http://www.general-data.com/hc)

*General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.*

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



Email: [talktous@general-data.com](mailto:talktous@general-data.com)  
 Web: [www.general-data.com/healthcare](http://www.general-data.com/healthcare)

General Data Healthcare, Inc.  
 4043 McMann Road, Cincinnati, OH 45245  
 844.643.1129

Service:  
 PH: 800.245.5029  
 Email: [customersupport@general-data.com](mailto:customersupport@general-data.com)