


SHURStain™ Manual Stainers



SS-3L
3-Position
Manual Stain Rack



SS-6L
6-Position
Manual Stain Rack




SS-12L
12-Position
Manual Stain Rack

Choose the Size that Best Suits Your Needs

These modular stainers are available by size to allow you to select the number of individual staining positions necessary to meet the needs of your lab or your individual or special staining application. They are produced in 3, 6, and 12 well configurations, and multiple units can be interlocked back-to-back to increase staining volume capabilities. Each SHURStain rack comes with a supply of documentation strips to record stain type, brand, lot#, etc. The number of data positions is proportional to the number of wells.

Lids Allow Easy Access with Minimal Solvent Loss

Sliding lids prevent evaporation by keeping stains covered.



Stainer Accessories

A digital timer attaches directly to the stainer with a convenient bracket. Green and white 250 ml Replacement Wells have sliding lids that open from either end to minimize evaporation.



Product Specifications:

Manual Stainer Rack	SS-3L	SS-6L	SS-12L
Dimensions: DxWxH (in)(mm)	(6x8x7)(160x203-x178)(in)(mm)	(6x15x6)(160x375x156)(in)(mm)	(6x28x6)(160x718x156)(in)(mm)

Catalog #	Description
SS-3L	SHURStain™ 3 Position Manual Stain Rack; includes 2 white 250ml wells (SS-WLW), 1 green solvent-resistant 250 ml well (SS-WLG), 1ea. 24-position plastic slide carrier (SS-SRG-24). Digital timer and a supply of label strips included (SS-LS-3L).
SS-6L	SHURStain™ 6 Position Manual Stain Rack; includes 4 white 250ml wells (SS-WLW), 2 green solvent-resistant 250 ml wells (SS-WLG), 24-position plastic slide carrier (SS-SRG-24). Digital timer and a supply of label strips included (SS-LS-6L).
SS-12L	SHURStain™ 12 Position Manual Stain Rack; includes 8 white 250ml wells (SS-WLW), 4 green solvent-resistant 250 ml wells (SS-WLG), 24-position plastic slide carrier (SS-SRG-24). Digital timer and a supply of label strips included (SS-LS-6L).
SS-SRG-24	Slide Rack, plastic, Grey, 24-position, solvent resistant
SS-LWL	Well, 250ml, w/ lid, white
SS-WLG	Well, 250ml. w/ lid, green, xylene resistant
SS-DT	Digital Timer, countdown, built-in memory, loud alarm, large display
SS-LS-3L	Label Strip, replacement strip for SS-3L with 3 data entry positions, 8", 12 ea/pk
SS-LS-6L	Label Strip, replacement strip for SS-6L with 6 data entry positions, 15", 12 ea/pk
SS-LS-12L	Label Strip, replacement strip for SS-12L with 12 data entry positions, 28", 12 ea/pk

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2018 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



Email: talktous@general-data.com
 Web: www.general-data.com/healthcare

General Data Healthcare, Inc.
 4043 McMann Road, Cincinnati, OH 45245
 844.643.1129

Service:
 PH: 800.245.5029
 Email: customersupport@general-data.com