

Histology Innovation for a NEW Generation

Consumables

SHURStor Slide & Cassette Storage File



**Proven Technology. Superior Reliability.
Protecting Every Life Story In Your Lab.**

- ✓ **Economical Storage For Your Cassettes and Slides**
SHURStor™ stores 20,000 slides or 5,000 cassettes/case
- ✓ **Space Saving Solution for Every Lab!**
SHURStor™ is a convenient portable lightweight file storage that can be stacked up to 10 files high.
- ✓ **Protects Slides & Cassettes for Long-term Storage**
SHURStor™ storage files are made of sturdy, rugged heavy-duty moisture-resistant cardboard that protects laboratory cassettes and slides for long-term storage of samples embedded in paraffin.

Product Specifications:

Dimensions: (in)	
SF-S	(11.50 x 11.50 x 3.50) (29.21 x 29.21 x 8.89)(in) (cm)
SF-B	(11.50 x 11.50 x 3.50) (29.21 x 29.21 x 8.89)(in) (cm)

How To Order:

SF-S	SHUR Stor ™ Slide Storage File w/4 drawers, capacity 2000 slides/file; 10 stackable files/cs.
SF-B	SHUR Stor ™ Tissue Block Storage w/8 drawers, capacity 500 cassettes/file; 10 stackable files/cs.

Related Products:

Consumables:	CLEAR View ™ Slide Holder, SHUR Mark ® One-Piece Cassettes, SHUR Track ™ Cassettes, SHUR Mark ® Microscope Slides, SHUR Mark ® Microscope Slide Folders / Boxes / Mailers
Instrumentation:	SHUR Mark ® Slide Labeler, SHUR Mark ® Cassette Labeler, SHUR Stain ™ H&E Slide Stainer, SHUR Stain ™ Manual Stainers, SHUR Dry ™ Slide Dryers, SHUR Clean ™ Fume Hoods


More Information

PH: 844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2014 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
Headquarters: 4354 Ferguson Drive, Cincinnati, OH 45245
TBS Facility: 3014 Croasdaile Drive, Durham, NC 27705

Sales Inquiries:
PH: 844.643.1129
Email: talktous@general-data.com
Web: www.general-data.com/hc

Service:
PH: 800.245.5029
Email: customersupport@general-data.com