

*Histology Innovation for a NEW Generation*

## Staining & Coverslipping

# Microscope Slide Drying Stand



**Proven Technology. Superior Reliability.  
Protecting Every Life Story In Your Lab.**

- ✓ **Durable and Chemical Resistant Stand**  
The polymer composition of the plastic resists staining and washes clean with ethanol, isopropanol or other organic solvents.
- ✓ **Space Saving Benchtop Design for Easy Addition to Microtomy Stations**  
Compact size is a great accessory to sit beside to SHURCut™ Microtome, TBS' **Flotation Work Stations**, SHURStain™ Manual or H&E Slide Stainers
- ✓ **Capacity To Accommodate Up To 77 Microscope Slides**
- ✓ **Angled Slots Permit Ease of Access to Individual Slides**
- ✓ **Allows for Rapid Viewing of Slide Labels**

# Product Specifications:

---

Dimensions: (in) (0.5 x 7 x 11.75) (in)

# How To Order:

---

SDS Slide Drying Stand, White, HDPE, stain resistant, 0.5" x 7" x 11.75"

# Related Products:

---

**Consumables:** CLEARView™ Slide Holder, SHURMark® Microscope Slides

**Instrumentation:** SHURStain™ Manual Stainers, SHURStain™ H&E Slide Stainer, Flotation Work Stations, SHURCut™ Microtomes


## More Information

**PH: 844.643.1129**

**[www.general-data.com/hc](http://www.general-data.com/hc)**

*General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.*

©2014 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.  
Headquarters: 4354 Ferguson Drive, Cincinnati, OH 45245  
TBS Facility: 3014 Croasdaile Drive, Durham, NC 27705

Sales Inquiries:  
PH: 844.643.1129  
Email: [talktous@general-data.com](mailto:talktous@general-data.com)  
Web: [www.general-data.com/hc](http://www.general-data.com/hc)

Service:  
PH: 800.245.5029  
Email: [customersupport@general-data.com](mailto:customersupport@general-data.com)