

SlideTrack™ ESPO

Direct to Slide Printer



■ On-Demand Slide Printing

The SlideTrack ESPO Slide Printer uses thermal transfer technology to print high quality 300 dpi text and barcodes directly on your slides at the microtomy station.

■ Automated Slide Printing That Meets Your Needs

The SlideTrack ESPO has a total capacity of 200 slides, with (2) 100-slide magazines; capable of printing 2 different types of slides. Three print modes allow your choice of single or batch printing. Additionally, the unique front loading feature lets you print worry-free on control slides!

■ Compact Footprint & Fast Speed

The small size allows for easy placement, and the SlideTrack ESPO prints slides at the rate of 3 seconds per slide (20 per minute).

■ Interface or Stand-Alone

Our experienced staff will assist with your installation and will ensure your lab is up and running smoothly.

Protect every life story™ in your lab.

Product Specifications:

Dimensions: WxDxH (in) [cm]	(4.7x9.4x9.8) [120x240x250]	Slide Size:	25mmx75mmx1mm, 26mmx76mmx1mm
Weight: (lbs) [kg]	(15.4) [7.0]	Frost Size:	15mm, 20mm, 25mm
Voltage:	100-240 VAC	Slide Capacity:	Two 100 slide magazines plus single slide front loading
Print Speed:	3 seconds per slide/20 slides per minute	Printer Interface:	USB (1), Ethernet, WiFi



Catalog #	Description	
G1-ESP-001W	SlideTrack™ ESPO Direct to Slide Printer, White <ul style="list-style-type: none"> • (1) Single-slide receiver • (1) 20-slide stacker • (2) 100-slide magazines • (1) ribbon (prints min. of 5000 slides) • (1) USB cable • Remote installation/label format support (G1-7DD-9000) 	
G1-ESP-S50	SlideTrack ESPO Slide Stacker: receives 50 printed slides	
G1-ESP-R025	Thermal Transfer Ribbon; 2 rolls/box, 5 boxes/cs, total of 10 rolls	
G1-ESP-MH100	SlideTrack ESPO Slide Magazine holder, holds 100 slides	

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com