

SlideTrack ESPO II Slide Printer



Fast, Powerful & Easy To Use Microscope Slide Printer

■ On-Demand Printing Directly On The Slide

Prints high quality 300 dpi text and barcodes directly on frosted or colored slides for positive specimen identification, minimized errors, and increased patient safety. The printed image is resistant to harsh histology lab chemicals.

■ Compact Footprint & Fast Speed

The small size allows for easy bench placement next to the microtome and water bath, and prints slides at the rate of 3 seconds per slide (20 per minute).

■ Dual Touch Screen Display

The color touch screen is configured for simultaneous display and operation modes offering intuitive controls.

■ Large Capacity With Single Slide Printing Capability

200 slide total capacity with two 100-slide magazines, additional single slide front loading feature allows for refill-free printing even with control slides when tissue is on the slide.

■ Built-In Barcode Scanner

An integrated barcode scanner on the front of the ESPO II saves bench space at the workstation by eliminating the need for a separate scanner.

■ Batch Or On-Demand Printing - Standalone or Connected To Your LIS

Supports either batch or on-demand slide printing, can operate connected to a single computer or our integration specialists can interface the ESPO II to your LIS.



Specifications

Part Number	ST-ESP-PTR-II
Dimensions	(W x D x H) 4.7" x 9.4" x 9.7" (119.4 mm x 238.8 mm x 246.4 mm)
Weight: lbs, kg	9.9 lbs (4.5 kg)
Print Resolution	300 dpi
Print Speed	3 seconds per slide / 20 per minute
Print Technology	Thermal transfer (requires thermal transfer ribbon)
Dot Density	12 dots/mm
Print Direction	0°, 90°, 180°, 270°
Solvent Resistance	Water / Alcohol / Xylene
Connectivity	USB, LAN, WLAN
Slide Capacity	Two 100-slide magazines + single slide front loading
Barcode Scanner	Integrated
Power Requirements	AC 90-240V

Accessories

ST-ESP2-MAG	ESPO II Slide Magazine - Holds 100 Slides
ST-ESP2-S50	ESPO II Front Loading 50 Slide Stacker
ST-ESP2-S20	ESPO II Front Loading 20 Slide Stacker
ST-ESP2-S1	ESPO II Manual Single Front Loading Tray
ST-ESP-R02-B	Thermal Transfer Ribbon with integrated spools; 1 box, 2 rolls/box
ST-ESP-R1-C	Thermal Transfer Ribbon with integrated spools; 1 case, 5 boxes/case, 2 rolls/box
ST-ESP2-RH	ESPO II Ribbon Cassette Holder Unit

Adhesion Microscope Slides – All Purpose & IHC



Matsunami APS
All Purpose Hydrophilic Slides



Matsunami TOMO®
IHC Hydrophilic Slides



Pre-Loaded Slide Cartridges

Get more info at www.general-data.com/slides

Value-Add Services



Installation & Set Up



Technical Support



**LIS Integration & Custom
App Development**



**Maintenance Plans &
Options**

For more information call us at 844-643-1129 or visit www.general-data.com/esp2