

# TempCheck<sup>o</sup>

Health Screening & Temperature Measuring Kit



## ■ Monitor the Health of Employees, Customers, & Visitors

Made by our partner SATO, the CT4-LX TEMPCHECK kit provides a ready-to-use solution for businesses to monitor the health of people entering their workplace, and safeguards the information within a **confidential logging system**.

## ■ Kit Contains Printer, Software, and Labels

The TEMPCHECK Kit is ready to go upon arrival, including printer, software, and 4"x2" labels. It also includes employee and customer/visitor versions of the TEMPCHECK software so the data is kept separately.

## ■ Pre-Loaded with COVID-19 Health Screening Questions

**HIPAA** compliant: entered health information is kept confidential within the logging system. As this is a stand-alone system, any temperature-reading device can be used.

## ■ CT4-LX Printer

With a 4.3-inch full-color touchscreen display and on-board intelligence, the CT4-LX seamlessly combines functionality, intelligence, and performance.

## ■ Equipped with TEMPCHECK Software

Easy-to-use, with prompts easily guiding the user to the next step. Includes health screening questions and temperature entry to monitor the health of employees and visitors, as well as a confidential logging system for information security.

## ■ Customization Available

If you want to scan to enter employee, patient, resident, or student information, or require other adjustments, just ask...customization is available for an additional fee.



# Product Specifications:

Software:	Employee & Vendor Versions of TEMPCHECK AEP software
Labels:	Direct thermal 4" x 2" labels, 4 rolls/cs
Hardware:	CT4-LX, direct thermal, 203 dpi
Catalog number:	CT4LXTEMPCHECK100



## CT4-LTx 4-inch Printer

So easy to use! With a 4.3-inch full-color touchscreen display and on-board intelligence, the CT4-LX seamlessly combines functionality, intelligence and performance.



## TEMPCHECK Software

Ready-to-use TEMPCHECK software includes COVID-19 health screening questions and temperature entry to monitor the health of employees and visitors, as well as a confidential logging system to safeguard the entered information.



## 4x2-Inch Labels

Kit includes 4 rolls of bright white, top-coated direct thermal labels with removable adhesive so labels can adhere to clothing. Labels are also able to fit into ID badge holders.

## More Information

844.643.1129

[talktous@general-data.com](mailto:talktous@general-data.com)

[web landing page address here](#)

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.  
4384 Ferguson Rd., Cincinnati, OH 45245  
844.643.1129

Email: [talktous@general-data.com](mailto:talktous@general-data.com)

Web: [www.general-data.com](http://www.general-data.com)

Service:

PH: 800.245.5029

Email: [customersupport@general-data.com](mailto:customersupport@general-data.com)